

BlueAQUA Prestained Protein Ladder

Catalog Number: PM019-0500N Size: 500 μ l
Catalog Number: PM019-0025N Size: 25 μ l
Lot Number: 15131037
Minimum Shelf Life: Mar-2023

Product Qualification

Specifications: The product must show 11 major protein bands when separated by SDS-PAGE electrophoresis on EVOGel with Tris-Glycine buffer, 12% Bis-Tris Gel with MES buffer, or 4-12% Bis-Tris Gel with MES buffer.

Results: Meet specifications.

Specification: Any pre-stained contaminant band must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.



5 μ l, EVOGel with Tris-Glycine-SDS buffer,



Yi Hsing, Yao

Yi-Hsing, Yao

Inspector, Quality

Issued on 2021/03/04