

Certificate of Analysis



0.25% Trypsin-EDTA, 1X

Catalog Number: CC508-0100

Lot Number: GL20151

Size: 100 mL

Minimum Shelf Life: Oct-2022

Storage: -30~-10 °C

Test	Specifications	Results
Osmolality	Check & Record	289
pH	Check & Record	7.4
Enzyme Activity	Check & Record	1268.67 U/ mL
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

A handwritten signature in black ink that reads 'Yi Hsing, Yao'.

Yi-Hsing, Yao
Inspector, Quality
Issued on 2020/10/30