

# Certificate of Analysis



## HEPES, 1M

Catalog Number: CC519-0100

Lot Number: GL19175

Size: 100 mL

Exp. Date: Aug-2021

Storage: 2 ~ 8 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	0.02 EU/mL
Osmolality	Check & Record	1744
pH	Check & Record	7.3
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.