



## 75T Cell Culture Flask

Cat. No.	Cell Growth Area(cm <sup>2</sup> )	Cap Style	Volume (mL)	Recommended Medium Volume(mL)	/Pack	/Case
PC273-0100	75	Sealed Cap	250	15-22.5	5	100
PC275-0100	75	Vent Cap	250	15-22.5	5	100

**Description:** Tissue Culture Treated Cell Culture Flasks, 250mL

**Sterile:** Yes

### Materials

Flask: PS (Polystyrene)

Cap: HDPE (Polyethylene)

Vent cap: 0.22µm hydrophobic filters

Simply cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures.

### Features

- Temperature Range: Approx.-86 to+64°C
- High clarity medical grade polystyrene material
- Sterilized by E-Beam SAL=10<sup>-6</sup>
- Non-pyrogenic, DNase/ RNase free, free, Non-Cytotoxicity
- Packaged in sterile, zip-sealable bags
- Two types of cap style : vent caps & plug seal cap
- Vent caps with 0.22µm hydrophobic filters for gas exchange without contamination
- Fostered writing and clear graduation
- Shelf Life: 3 Years after date of production (ensure package is in good)
- Clear lot number for batch traceability
- Stackable

