

# Certificate of Analysis



## Ham's F-12

Catalog Number: CC119-0500

Lot Number: GL19534

Size: 500 mL

Exp. Date: Mar-2021

Storage: 2 ~ 8 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.007 EU/mL
Osmolality	Check & Record	284
pH	Check & Record	7.8
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.