

## UltraTRAX Transfection Reagent

Catalog Number.MB508-1000 Size: 1000  $\mu$ L,

Catalog Number.MB508-0100 Size: 100  $\mu$ L,

Lot Number: 95156100

Store at -20°C


### Product Qualification

Appearance: Clear, colorless liquid.

Cell: 293T ( $5 \times 10^5$  cell in 6 cm dish)

Transfection condition:

FLAG-mCherry plasmid 2.5  $\mu$ g with UltraTRAX Transfection Reagent 7.5  $\mu$ l



All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.