

Certificate of Analysis



MBead Bacteria Genomic DNA Kit

Catalog Number: NA010-0004

Size: 4 reactions

Lot Number: 80110300

Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	Lysis Buffer	Wash Buffer	Release Buffer	Magnetic Bead
pH	4.5	5.0	8.5	---
accredited	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	DH5α	300µl	30.8	√	√
2	DH5α	300µl	31.4	√	√
3	DH5α	300µl	32.6	√	√
4	DH5α	300µl	31.3	√	√
5	DH5α	300µl	32.5	√	√
6	DH5α	300µl	30.7	√	√

Result: Meets specification.

Storage condition: MBead Bacteria Genomic DNA Kit should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses