

## PiNK Plus Prestained Protein Ladder

Catalog Number. PM005-0500

Size: 500  $\mu$ l

Catalog Number. PM005-0025

Size: 25  $\mu$ l

Lot Number: 09915021

### Product Qualification

Specifications: The product must show 11 protein bands when separated by SDS PAGE electrophoresis on an OnePAGE Sonic Gel with Tris-Glycine-SDS buffer, 12% Bis-Tris Gel with MOPS buffer, and 4-12% Bis-Tris Gel with MES buffer.

Results: Meet specifications.

Specification: Any contaminant band must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.



5  $\mu$ l, OnePAGE Sonic Gel with Tris-Glycine-SDS buffer

*Yi Hsing, Yao*

---

Yi-Hsing, Yao

Inspector, Quality

Issued on 2020/12/22