

Certificate of Analysis



PiNK Plus Prestained Ladder

Catalog Number: PM005-0500

Lot Number: 00145027

Size: 500 µl

Minimum Shelf Life: May-2022

Product Qualification

Specifications: The product must show 11 bands when separated by SDS PAGE electrophoresis on an OnePAGE Sonic Gel with Tris–Glycine buffer, 12% Bis–Tris Gel with MOPS buffer, or 4-12% Bis–Tris Gel with MES buffer.

Results: Meet specifications.

Specification: Any contaminant band must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.



5ul, OnePAGE Sonic Gel with Tris/ Glycine/ SDS buffer

Yi-Hsing, Yao

Inspector, Quality

Issued on 2020/05/19