

BlueRAY Prestained Protein Ladder

Catalog Number: PM006-0500N

Lot Number: 80106802

Size: 500 μ l

Product Qualification

Specifications: The product must show 10 major protein bands when separated by SDS-PAGE electrophoresis on a EVOGel with Tris-Glycine-SDS buffer, Tris-glycine-SDS PAGE, 12% Bis-Tris Gel with MOPS buffer, or 4-12% Bis-Tris Gel with MES buffer.

Results: Meet specifications.

Specification: Any pre-stained contaminant band must be less than 10% in intensity to the nearest major band.


Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



5 μ l, EVOGel with Tris-Glycine-SDS running buffer