

## BLUeye Prestained Protein Ladder

Catalog Number: SP007-0600

Lot Number: 01061131

Size: 300  $\mu$ L X 2 vials

### Product Qualification

Specifications: The product must show 12 major protein bands when separated by SDS-PAGE electrophoresis on a EVOGel with Tris-Glycine-SDS running buffer.

Results: Meet specifications.

Specification: Any pre-stained contaminant band with molecular weight below 245 kDa must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



5 $\mu$ l, EVOGel with Tris-Glycine-SDS running buffer

*Yi Hsing, Yao*

---

Yi-Hsing, Yao  
Inspector, Quality  
Issued on 2020/08/07