

Certificate of Analysis



OnePCR

Part No./Catalog No. MB203-0004

Lot Number: MB20120025

Minimum shelf life: Feb-2022

Size: 4 Reaction

Product Qualification

Specifications: The product must show the DNA amplification by PCR reaction. The bands when separated by 1.0% agarose electrophoresis with 1X TBE Buffer

Results: Meet specifications.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



The DNA was amplified from p35S-GFP plasmid.

Gel was explored under the Blue light.

Handwritten signature of Hsien-Yi Lin.

Hsien-Yi, Lin

Inspector, Quality

Issued on 2021, Feb